



SOMNODENT PRODUCT CHART

| | PREMIUM LINE | SIGNATURE LINE | | | STANDARD LINE | | |
|----------------------------------|--|---|---|---|---|---|--|
| | | | | | | | |
| | AVANT | FUSION | FLEX | CLASSIC | AIR / AIR+ | HERBST ADVANCE | MORNING REPOSITIONER |
| WARRANTY | 3 YEARS | 3 YEARS | 3 YEARS | 3 YEARS | 1 YEAR | 2 YEARS | 1 YEAR |
| PRODUCT FEATURES | Slim and durable, the Avant is the only milled oral device featuring their soft, patented Bflex liner for optimal patient comfort. Their unique titration system makes calibration easier than ever. | An innovation on the original design with the fusion of interchangeable wings and screw calibration, offering quick advancement at 0.1mm increments and an 8.5mm range. | Features a soft inner liner made with a unique polymer material, SMH BFlex, which maintains retention and offers premium patient comfort. | The original design in customized oral devices, the Classic is held in place using medical grade stainless steel ball clasps. | Comfort and high-quality features in an affordable device, optimizes low reimbursement scenarios. Available with ball clasps or standard Poly+ liner for retention. | Visual calibration indicator allows greater control over therapy. Less resets needed due to 8 mm range of calibration. | Hands-free bite deprogrammer designed to return the patient's mandible back to pre-treatment centric position. |
| RECOMMENDED CLEANING AND STORAGE | Rinse with cold or lukewarm water, soak with a SomTab in a container of water, rinse, then dry and store in a dry container with the lid off. | Rinse with cold or lukewarm water, soak with a SomTab in a container of water, rinse, then store in a container of clean water. | Rinse with cold or lukewarm water, soak with a SomTab in a container of water, rinse, then store in a container of clean water. | Rinse with cold or lukewarm water, soak with a SomTab in a container of water, rinse, then store in a container of clean water. | Rinse with cold or lukewarm water, soak with a SomTab in a container of water, rinse, with lukewarm water, dry and store in a dry container. | Rinse with cold or lukewarm water, soak with a SomTab in a container of water, rinse, then store in a container of clean water. | Rinse with cold or lukewarm water, dry and store dry. |
| TURNAROUND TIME | 14 WORKING DAYS | 14 WORKING DAYS | 14 WORKING DAYS | 14 WORKING DAYS | 14 WORKING DAYS | 14 WORKING DAYS | 14 WORKING DAYS |

SOMNOMED ORAL DEVICE RECOMMENDATIONS

✔ Device is specially made for this patient type
 ✔ Device will work very well for this patient type
 ✔ Device will work for this patient type

| PATIENT CONDITION | SUGGESTED FUNCTION | PREMIUM LINE | FUSION | FLEX | CLASSIC | AIR / AIR+ | HERBST ADVANCE | MORNING REPOSITIONER |
|--|---|--------------|--------|------|---------|------------|----------------|----------------------|
| Open and closing mouth is a priority for the patient | Allows mouth to open / close more freely without metal hinges | ✔ | ✔ | ✔ | | ✔ | | ✔ |
| Larger Tongue | Increased tongue space for patient | ✔ | | | | | ✔ | ✔ |
| Bruxism / TMD Sufferer | Increased tongue space for patient | | | | | | ✔ | ✔ |
| Smaller Mouth | Shorter adjustment arms for compatibility / Fusion is slimmer with less bulk due to prefabricated wings | ✔ | ✔ | | | | ✔ | ✔ |
| Smaller Mouth & Bruxism/TMD Sufferer | Shorter adjustment arms for compatibility and durability for clenching & grinding | | | | | | ✔ | ✔ |
| Mouth Breathers | Option of anterior opening for capability of mouth breathing while device is worn | ✔ | ✔ | | ✔ | | ✔ | ✔ |
| Extreme vertical opening / Patient with large vertical | Option of elastic retention (ER) hooks | | ✔ | ✔ | | | ✔ | ✔ |
| Transportation patient / Positional Therapy Patient | DentiTrac not only allows for compliance tracking, but sleep position as well | | ✔ | ✔ | | | ✔ | ✔ |

